



i-Fix Implant

“i-Fix implants: Made in India for the World”

Dr Mohammad Shahnawaz

Kamal Medtech Who we are?

Global Footprint



Unit 1 - Manufacturing Facilities at Yamunanagar -200Km north of Delhi



*Unit 2 - Casting & Machining Facilities
at Yamunanagar*



Unit 3 - New Facility at IMT Faridabad, NCR Delhi



Kamal Medical Devices Brands



Stenoflex
Sirolimus Eluting Coronary Stent System

Everoshine
Everolimus Eluting Coronary Stent System

i-FIXTM
— POWERING SMILES —

The Beginning.....

Research Partners



Indian Institute of Technology, Delhi

iFix Dental Implant has been Conceptualized, Designed & Developed by Indian Institute of Technology (IIT), Delhi in collaboration with Maulana Azad Institute of Dental Sciences (MAIDS), Delhi



Maulana Azad Institute of Dental Sciences
New Delhi

iFix Dental Implant has been Clinically studied at Maulana Azad Institute of Dental Sciences, New Delhi



Council of Scientific and Industrial Research

Under the Patronage of Council of Scientific & Industrial Research (CSIR) - NMITLI Program



Manufacturers of high quality medical devices of global standards

Certifications: MDR 2017 & ISO 13485



FORM MD-9

[See sub-rule (1) rule 25]

Licence to Manufacture for Sale or for Distribution of Class C or Class D medical device

Licence Number: MFG/MD/2019/000057

- M/s KAMAL ENCON INDUSTRIES LIMITED, 56 INDUSTRIAL ESTATEYAMUNA NAGAR, -, Haryana (India) - 136001 Telephone No.: 1123351365 FAX: 1123721657 has been licenced to manufacture for sale or for distribution the below listed medical device(s) at the premises situated at M/s KAMAL ENCON INDUSTRIES LIMITED, Plot No. 917, IMT Faridabad, Sector- 68, Faridabad, Haryana (India) - 121001 Telephone No.: 11-23723158, 11-23351369 FAX: 11-23721657
- Details of medical device(s) [Annexed]
- The names, qualifications and experience of the competent technical staff responsible for the manufacture and testing of the above mentioned medical device(s); As per records maintain by the manufacturer
- This licence is subject to the provisions of the Medical Devices Rules, 2017 and conditions prescribed therein.

ANNEXURE

S.No.	Details Of Device(s)
1	<p>Generic Name:Dental Implants with Abutments & Prosthesis Model No.:NIL Intended Use:Dental Implants are intended for surgical placement in the upper or lower jaw to provide a means for prosthetic attachment in single tooth restorations and in partially or fully edentulous spans with multiple single teeth utilizing delayed loading, or as a terminal or intermediary abutment for fixed or removable bridgework, and to retain overdentures. Class of medical device:Class C Material of construction:Titanium Alloy, Ti-6Al-4V Dimension(if any):Diameter (mm): 3.00, 3.50, 4.00, 4.50, 5.00 & 5.50 AND Length (mm): 7.50, 8.50, 10.00, 11.50, 13.00 & 15.00 Shelflife:60 months Sterile or Non sterile:STERILIZED Brand Name:(if registered under the Trade Marks Act, 1999):iFix, MultiFix, SFix, CFix</p>

Place:
Date 29-May-19

S ESWARA
REDDY
Central Licensing Authority

| evaluate | assess | certify |



assessors of quality

Certificate No. MD-QMS/91/R/2609

This is to certify that

Kamal Encon Industries Ltd.

Plot No. 917, Sector 68, IMT Faridabad – 121 001, Haryana, India

has been found to conform to the requirements of
Medical Devices - Quality Management System Standard

ISO 13485:2016

This certificate is valid for the following scope:

Design, Development, Manufacture and Sales of:

- Interventional Cardiology Products:**
 - Sirolimus Eluting Coronary Stent System
 - Everolimus Eluting Coronary Stent System
 - Dental Implants and Prosthesis**
- Manufacture and Sales of Surgical Kits**

Initial Certification : 27/07/2019
Valid until : 26/07/2022



Authorized Signatory

The authenticity of this document may be verified by visiting www.zenithqa.com. Lack of fulfillment of conditions as set out in the Certificate Agreement may render this certificate invalid, any alteration, forging or falsification of the certificate or production of the same for purposes not intended by the holder shall be liable.

Accreditation Body: National Accreditation Board for Certification Bodies (NABCB)
Quality Council of India (QCI), 2nd Floor, 1st Main Road, Highways Building, Sector 28, Gurgaon,
New Delhi - 110 002, India
www.nabcb.org.in

Issuing Authority: Zenith Quality Assessors Pvt. Ltd.
308, 4th Floor, Sector 46, Gurgaon, Haryana - 122 002, India
www.zenithqa.com



CERTIFICATE

EC Certificate No. 1434-MDD-344/2021
Full Quality Assurance System
Directive 93/42/EEC concerning medical devices

Polish Centre for Testing and Certification certifies
that the quality assurance system in the organization:

Kamal Encon Industries Limited
Plot No. 917, Sector-68,IMT, Faridabad 121001,
Haryana, INDIA

for the design, manufacture and final inspection of
medical devices, class IIb

Dental Implants & Prosthesis

The list of medical devices covered by this certificate is provided in the Annexes 1-13

complies with requirements
of Annex II (excluding Section 4) to Directive 93/42/EEC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 25.05.2021 to 27.05.2024

The date of issue of the Certificate: 25.05.2021

The date of the first issue of the Certificate: 25.05.2021




Issued under the Contract No. MD-80/2020
Application No: 138/2020
Certificate bears the qualified signature.
Warsaw, 25/05/2021
Module H2/3/4/5

Vice-President


IIT Delhi Conceptual Design - led to Indian and US Patent




US Patent No. 9,833,300




Assembly



Implant



Abutment



Screw

GRADE 5 ELI



Properties		Advantages
Tensile strength	860 MPa minimum	Extremely high strength/weight ratio
Elongation	10%	Improved ductility and fracture toughness under static and dynamic loads
Oxygen content	0.13% maximum	
Biocompatibility	Excellent	Direct structural and functional connection between bone and implant surface


Ti Grade 5 ELI




- Its properties of biocompatibility, corrosion resistance, high strength, and toughness make it especially suitable for implants and other medical devices.
- ELI is also referred to as Grade 23 Titanium and ASTM F136.
- Ti6Al4V ELI contains reduced levels of oxygen, nitrogen, carbon and iron. ELI is short for “Extra Low Interstitials”, and these lower interstitials provide improved ductility and better fracture toughness for the Ti6Al4V ELI material.

Properties		Advantages
Tensile strength	860 MPa minimum	Extremely high strength/weight ratio
Elongation	10%	Improved ductility and fracture toughness under static and dynamic loads
Oxygen content	0.13% maximum	
Biocompatibility	Excellent	Direct structural and functional connection between bone and implant surface

The Design Premise






FE Meshing



Bottom all degree of constrains

Boundary Condition

Comparison of different thread types under compression, Tension						
S. No.	PARAMETER & LOCATION	Classical V	Reverse Buttress	Buttress	Buttress & Micro threads	Reverse Buttress & Micro threads
I) COMPRESSION		Units-Mpa				
A)	Implant Abutment interface	41.49	36.51	30.94	38.49	47.08
B)	Crestal bone level	29.71	11.94	20.46	18.48	23.73
C)	Cortico-cancellous interface	30.44	24.39	27.25	18.85	34.88
I) TENSION						
A)	Implant Abutment interface	40.47	37.57	33.95	38.44	47.08
B)	Crestal bone level	29.73	44.38	20.44	18.48	23.74
C)	Cortico-cancellous interface	30.99	63.34	27.25	22.29	34.88
III) MOMENT						
A)	Implant Abutment interface	425.34	421.62	278.94	396.09	339.45
B)	Crestal bone level	227.03	120.68	139.97	70.27	49.55
C)	Cortico-cancellous interface	9.59	9.18	30.39	8.57	11.23



Design of CSIR-NMITLI Dental Implant System




Implant

Abutment

Screw

Ball Abutment


Impression Coping



Implant Abutment Assembly




Implant Ball Abutment Assembly



Implant Healing Screw Assembly


Analysis and Testing as per ISO 14801 & 13498



Loading direction


Table 3: Fatigue test results

Serial No.	Load value (in percent of the static load to failure)	Life (No. of cycles)
1.	80 %	3,838
2.	70 %	18,943
3.	60 %	36,726
4.	50 %	89,840
5.	40 %	5,000,000




300 N


350 N



400 N




Fatigue Testing ISO 14801



Torsion Testing ISO13498


A comparative study with commercially available Implants




- **Animal Trials on rabbit model**
- **Physical Testing at IIT D**
- **Micro CT Analysis**
- **Histomorphometric Study**
- **Sterility Testing**
- **Stability testing**



Dental Implant Designing Evolution & Validation of Designs by FEM



Micro Thread
Improved biomechanical performance at Buttressbone-implant Thread interface.



Elastic Properties Ascribed to Materials Used In The Models

Materials	Young Modules (GPa)	Poisson's ratio
Cortical bone	15	0.3
Cancellous bone	1.5	0.3
Titanium	110	0.35





Fig. 3 Comparing different thread types under moment

- DYNAMIC FATIGUE TEST:
DONE AT IIT,DELHI ON
CETR-UMT-2, BRUKER
- Implant abutment
assembly with load
applied 30° to vertical for
fatigue testing (as per FDA
norms)
- • 2 MILLION CYCLES - 1
SAMPLE
- • 5N -800 N
COMPRESSION –
COMPRESSION TEST



Shri Ram Institute of Industrial Research



In Vivo Implantation study of specially prepared Micro implants were done on femur of White New Zealand rabbits and evaluated after 6 and 12 weeks. After necropsy histological, micro CT and Histomorphometric evaluation was done which revealed no adverse effects and an extremely favourable bone implant contact.



MICRO CT ANALYSIS




**Shree Chitra Tirunal Institute for Medical sciences
Technology (Bio Medical Technology Wing),
Thiruvananthapuram**

Sample Holder in Position for Analysis




RESULTS OF MICRO CT IMAGING FOR SAMPLE ON DAY 84/ 12 WEEKS (SAMPLE MNC3)




Medullary space

Trabecular projections





Prominent trabecular projections seen as strands extending from the inner aspect of the surrounding cortical bone towards the implant surface engaging it. High concentration of trabeculae seen near the area where the neck of the implant engages the cortical bone.

6 weeks



12 weeks



S.NO	SPECMEN/ WEEKS	MEAN DISTANCE FROM CRESTS IN RECTANGLE	PERCENTAGE BONE FORMATION
1	12 weeks	415.62 µm	80.5%
2	6 weeks	424.63µm.	25.6%



SHRIRAM INSTITUTE FOR INDUSTRIAL RESEARCH

(A unit of Shriram Scientific and Industrial Research Foundation)

An ISO - 9001:2008 Certified Institute

TEST CERTIFICATE

000358830

Issued to :
INDIAN INSTITUTE OF TECHNOLOGY, DELHI
DEPT. OF MECHANICAL ENGINEERING
HAUZ KHAS
NEW DELHI - 110016

J.O.No. 502-111-1839
Reg.No. 1370083
Date 21-02-2015
GC-01 (REV-04)
Your Ref.No. -



Kind Attn: DR NARESH BHATNAGAR , PROFESSOR & ASSO.
DEAN [R&D]

Sample Particulars :

One sample of Dental Implant labelled as Lot No. 1 DT 17.01.2015 was received.

"The sampling was not carried out by Shriram Institute for Industrial Research. The sample details provided in test certificate are based on declaration by the party."

TEST RESULTS

S.No	TEST	RESULT	PROTOCOL
1.	Sterility test	Sample complies with the test for sterility	As per guidelines of USP-37


NOTE: Party asked for above test only.

DOR: 06-02-15
DOC: 20-02-15




Indian & US Patent






**INTELLECTUAL
PROPERTY INDIA**
PATENTS | DESIGNS | TRADE MARKS
GEOGRAPHICAL INDICATIONS



भारत सरकार
GOVERNMENT OF INDIA
पेटेंट कार्यालय
THE PATENT OFFICE
पेटेंट प्रमाणपत्र
PATENT CERTIFICATE
(Rule 74 Of The Patents Rules)

क्रमांक : 011117822
SL No :



पेटेंट सं. / Patent No. : 324068


आवेदन सं. / Application No. : 2243/DEL/2007


फाइल करने की तारीख / Date of Filing : 26/10/2007

पेटेंटी / Patentee : COUNCIL OF SCIENTIFIC & INDUSTRIAL RESEARCH

प्रमाणित किया जाता है कि पेटेंटी को उपरोक्त आवेदन में यथाप्रकटित DENTAL IMPLANT SYSTEM नामक आविष्कार के लिए, पेटेंट अधिनियम, 1970 के उपबंधों के अनुसार आज तारीख 26th day of October 2007 से बीस वर्ष की अवधि के लिए पेटेंट अनुदान किया गया है।

It is hereby certified that a patent has been granted to the patentee for an invention entitled DENTAL IMPLANT SYSTEM as disclosed in the above mentioned application for the term of 20 years from the 26th day of October 2007 in accordance with the provisions of the Patents Act, 1970.





अनुदान की तिथि : 31/10/2018
Date of Grant :

पेटेंट नियंत्रक
Controller of Patent

Signature

टिप्पणी - इस पेटेंट के नवीकरण में गिरावट, यदि इसे नवापन गलत नाम है, 26th day of October 2009 को जी। जगन्ने शंकर प्रसेक वर्ष में जारी दिन रच होगी।
Note - The fees for renewal of this patent, if it is to be maintained will fall / has fallen due on 26th day of October 2009 and on the same day in every year thereafter.



(12) United States Patent Verma et al.

(10) Patent No.: **US 9,833,300 B2**
(45) Date of Patent: **Dec. 5, 2017**

- (54) **DENTAL IMPLANT SYSTEM**
- (75) Inventors: **Mahesh Verma**, New Delhi (IN);
Naresh Bhatnagar, New Delhi (IN);
Abhinav Sood, New Delhi (IN);
Farukh Faraz, New Delhi (IN);
Kshitij Sharma, New Delhi (IN);
Gedela V. Rao, New Delhi (IN); **Palani S. Kumar**, New Delhi (IN); **Shankar Iyer**, Elizabeth, NJ (US)

(58) Field of Classification Search
CPC A61C 8/0048- 8/0078
(Continued)

(56) **References Cited**
U.S. PATENT DOCUMENTS
5,588,838 A * 12/1996 Hansson A61C 8/0022
433/173
6,149,432 A * 11/2000 Shaw et al. 433/174
(Continued)

- (73) Assignee: **Council of Scientific & Industrial Research**, New Delhi (IN)
- (*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 188 days.

FOREIGN PATENT DOCUMENTS
WO WO 03/030767 A 4/2003
WO WO 2006/081239 A 8/2006
WO WO 2007/022655 A 3/2007

- (21) Appl. No.: **12/739,654**
- (22) PCT Filed: **Oct. 23, 2008**
- (86) PCT No.: **PCT/IN2008/000700**
§ 371 (c)(1),
(2), (4) Date: **Jan. 28, 2011**
- (87) PCT Pub. No.: **WO2009/054005**
PCT Pub. Date: **Apr. 30, 2009**

OTHER PUBLICATIONS
Eskitascioglu et al., "The influence of occlusal loading location on stresses transferred to implant-supported prostheses and supporting bone: a three-dimensional finite element study," Feb. 2004, Journal of Prosthetic Dentistry, vol. 91, pp. 144-150.*

- (65) **Prior Publication Data**
US 2011/0117522 A1 May 19, 2011
- (30) **Foreign Application Priority Data**
Oct. 26, 2007 (IN) 2243/DEL/2007

Primary Examiner — Matthew Nelson
(74) Attorney, Agent, or Firm — Locke Lord LLP; Daniel J. Fiorello; Joshua L. Jones

- (51) Int. Cl.
A61C 8/00 (2006.01)
- (52) U.S. Cl.
CPC **A61C 8/005** (2013.01); **A61C 8/0022** (2013.01); **A61C 8/0025** (2013.01);
(Continued)

(57) **ABSTRACT**
A screw type dental implant system (1), having a dental implant fixture (2), a multifunctional component (3) and an abutment screw (4). The screw shape dental implant fixture has an external surface having buttress threads on the body and micro threads at the collar. This combination provides the advantages of: improved biomechanics at the implant abutment interface, self-tapping nature to the implant, and minimizing the stresses at the crest of the bone leading to decrease resorption of crestal bone. The multifunctional component serves the purposes of implant mount, impression analog and final abutment and has a single prosthetic platform so one component is compatible with different implant dimensions, which minimizes the inventory needed for the implant system and allows easy handling of the
(Continued)

Kamal Medtech Facility

Facility Area - 50000 Sqft
State of the Art Clean Rooms
2017 & ISO 13485
High-tech QC Labs
Equipments
Well Trained Man Power

Controlled Area – 20000 Sqft
Compliant with MDR

World Class Mfg

Best-in-Class Technologies



Key Features



- Single Prosthetic Platform for all Diameters **hence inventory is minimized**
- Having the benefit of Platform Switch in all connections **hence better soft tissue aesthetics**
- Conical connection 12 degree morse taper which prevents **micro leakage**
- Anatomical tapered design
- Rough Surface – **for better hard and soft tissue integration**
- SLA Surface for **better osseointegration**
- Smart Surgical Kit – **Comprehensive, Compact and easy to use and caters to most surgical & prosthetic requirement.**

Optimized thread designs




- Buttress thread engineered for superior Primary stability
- Coronal bevel for better soft tissue integration and inbuilt platform shift.
- Reduction of stresses in bone at both Crestal & cancellous level

Topography evaluation by SEM of marketed dental implants




- Large grit sand blasted
- 3 Stage acid etching process
- Roughness is 1.2-2.0 Ra value on Buttress Thread
- Extensive multi-stage cleaning removes undesired residues ,yielding a contamination-free surface.

Magnification 2500




Magnification 5000



Implant Connection - Conical

- Internal Hex design-Rigid connection
- Stable biological seal
- 6 prosthetic positions
- 12° Morse taper-Creates an optimal seal against bacterial leakage
- Precise fit
- 2.25mm Hex connection diameter



IMPLANT DIMENSIONS

Implant Dimensions & Drilling Protocol

Size Table

Length (mm)

Dimensions	6.5	7.5	8.5	10	11.5	13	15
3.5	×	×	✓	✓	✓	✓	✓
4	✓	✓	✓	✓	✓	✓	✓
4.5	✓	✓	✓	✓	✓	✓	✓
5	✓	✓	✓	✓	✓	✓	×

Diameter (mm)



3.5 mm



4 mm



4.5 mm




5 mm

SURGICAL KIT



SURGICAL KIT



Amorphous DLC coated drills for longer shelf life. Length specific laser marking and stoppers for better control while drilling

Drilling Sequence as per speed*



Lance Drill
1200 rpm



Lindeman Drill
1000 rpm



2.0 drill
1000 rpm



Tapered Drill
600-800 rpm



Cortical Drill
600-800 rpm




Implant Driver
30-50 Ncm

*WITH COPIOUS IRRIGATION

i-Fix Pro

The unique implant design for high primary stability and better esthetic outcome



Double lead threads The expanding taper with double lead thread design of i-Fix PRO offers high initial stability, a gradual bonecondensing property, and adjustable implant orientation for optimal prosthetic connection angulation. In detail, the 1.3mm spacing of the double start threads advances the implant 2.6mm with each rotation. In contrast to i-Fix implants advance approximately 0.9mm with each rotation.

The coronal region of the i-Fix PRO 4/4.5&5mm diameter implants are backtapered, a design that allows for maximum bonevolume around the implant platform.

Implant body design The expanding tapered body is designed to attain high stability even in compromised bone situations.

Reverse-cutting flutes Running along i-Fix PRO are two reverse-cutting flutes (length of the flutes depends on implant length). Insertion of the implant counterclockwiseengages the cutting capability of these flutes, thereby gradually widening the osteotomy and allowing for continued insertion.

Apical drilling blades facilitate insertion into minimally prepared sites.designThe apical region of i-Fix Pro features twin cutting blades that



PROSTHETIC PORTFOLIO



TEMPORARY ABUTMENTS



- Available in Hex & Non-Hex
- Retentive grooves to engage acrylic temporary materials.
- For Single unit (Hex) and Multiple unit (Non-Hex) temporizations.
- For screw-retained single- and multiple-unit restorations.
- Extended length designed for customization per clinical requirements



Hex



Non-Hex

Temporary Abutment Hex	TAH4001
Temporary Abutment Non-Hex	TANH4001

TEMPORARY ABUTMENTS (IMPLANT LEVEL)



Place the abutments



Fabricate temporary using vacuum formed stent



Drill holes in the stent and remove material to access the screws



Remove provisional & finish it




Screw the bridge and seal the access holes with composite

Temporary Abutment
Titanium



HEALING ABUTMENTS



Diameter	Collar	Length	Article code
4.15	1	2	HA4012
4.15	1	4	HA4014
4.15	2	3	HA4023
4.15	2	5	HA4025
4.15	2	6	HA4026
4.15	3	3	HA4033
4.15	3	5	HA4035
4.15	3	6	HA4036
5.2	1.5	3	HA5153
5.2	1.5	4	HA5154
5.2	2.5	3	HA5253
5.2	3.5	3	HA5353

Healing Abutment :

- Use Hex driver
- Tightening torque : 5-8Ncm

IMPRESSION COPING – OPEN/ CLOSED



Open Tray Impression Coping :

- Custom tray is used for pick up impression taking.
- Use 1.2 Hex driver
- Packing contents: Impression Coping + Guide Pin

Collar	4.15	5.2
1mm	OTI4001	-
1.5mm	-	OTI5001
2mm	OTI4002	-
2.5mm	-	OTI5002



Closed Tray Impression Coping

- Custom tray is used for pick up impression taking.
- Use 1.2 Hex driver
- Packing contents: Impression Coping + Guide Pin
- Available as Long and Short.

Product Description	Article Code
Closed Tray Impression Coping- Long	CTICL
Closed Tray Impression Coping- Short	CTICS
Closed Tray Impression Screw- Long	CTISL
Closed Tray Impression Screw- Short	CTISS

SHOULDER ABUTMENT



Cement Retained | Single or Multiple

Margin Height Selected
based on Soft Tissue
Depth

Benefits

- Packaging unit: Abutment+ Abutment screw
- Tightening Torque: 35Ncm
- Available in multiple collar heights and 2 diameters.
- Available as Hex & Non-Hex.
- Versatile features make this abutment indicated in multiple clinical scenario.



4.15 Hex



4.15 Non-Hex



5.2 Hex



5.2 Non-Hex


Product Description	Hex	Non-Hex
S Abutment 4.15 x1	SAH4001	SANH4001
S Abutment 4.15 x2	SAH4002	SANH4002
S Abutment 4.15 x3	SAH4003	SANH4003
S Abutment 5.2 x1.5	SAH5001	SANH5001
S Abutment 5.2 x2.5	SAH5002	SANH5002
S Abutment 5.2 x3.5	SAH5003	SANH5003

ANGLED ABUTMENT




Pre-fabricated Titanium abutment for single retained restorations

- Margin designed to profile natural soft tissue contours
- 15 & 25 degrees angulation
- Indicated in anterior esthetic zone
- Comprehensive selection of varying margins and angulations minimizes chair-side adjustments.
- Use Hex Driver.
- Packaging contents: Abutment + Abutment Screw.
- Tightening torque: 35Ncm
- Length of abutment is 12.84mm



Product Description	15°	25°
Angled Abutment 4.15 x 1	15A4001	25A4001
Angled Abutment 4.15 x 2	15A4002	25A4002
Angled Abutment 5.2 x 1.5	15A5001	-
Angled Abutment 5.2 x 2.5	15A5002	-

UCLA ABUTMENT



UCLA with Hex 4.15x1mm	UCH4001
UCLA with Hex 5.2x1.5mm	UCH5001
UCLA Nonhex 4.15x1 mm	UCNH4001
UCLA Nonhex 5.2x1.5	UCNH5001
Abutment screw	AS01

Versatile:

- Custom abutment for all indications using the lost wax casting technique which can be used for both cement or screw retained restorations.
- For screw-retained restorations directly to the implant
- Burn-out cylinder for wax application
- For single (hex) or multiple (non-hex) units.
- Used to correct implant angulation and adaptable to limited inter-occlusal dimensions
- Custom abutment for the small labs without the purchase of capital equipment.
- Length of abutment is 11mm (8.4 mm from collar)

CUSTOMIZED SOLUTION : TIBASE ABUTMENT



Ti-Base Hex



Ti-Base Non-Hex

Ti Base Abutment:

- Titanium bonding base for the production of customized abutments made of zirconia for i-Fix Implants.
- No requirement of sandblasting
- CAD/CAM Zirconia superstructure bonded to titanium base
- Indicated for anterior esthetic zone.
- Screw-retained single tooth restorations: CAD/CAM zirconium dioxide superstructure and directly veneered crown bonded to titanium base
- Cemented single/ bridge tooth restorations: CAD/CAM zirconium dioxide superstructure bonded to titanium base, with full-ceramic crown
- Packaging unit: Tibase abutment + 2 abutment screws
- Tightening Torque: 35 Ncm

Ti Base abutment Hex	TBH40
Ti Base abutment Non-Hex	TBNH40

Will be available for Digital
Prosthetic workflow post April
2022


CUSTOMIZED SOLUTION :PREMILL BAR



- These are used to fabricate customized abutments
- Available in 2 diameters 10 & 14mm & total length of 28mm
- Intact internal connection.

Pre-milled bar 10mm	PB10
Pre-milled bar 14mm	PB14

Digital Workflow Solution : HScan Body & Digital Analog



Digital Analog



H scan Body


H-Scan Body :

- Use to take Digital Impressions & enables the digital workflow to be more efficient.
- Available in 2 diameters 4.15 & 5.2mm
- Available in 3 lengths 4mm, 6mm & 8mm
- Titanium-made – autoclavable and visible in X-ray scans

Full Arch Solutions: Over Denture Solution



Ball Abutment



Versatile:

- Packing unit: Ball Abutment + Abutment Holder
- Ball Abutment Cap package contains: Ball Cap, Spacer and O Rings
- Tightening torque: 25Ncm
- Ball cap insertion and removal tool is available for changing o rings

Ball Abutment 1mm collar	BA01
Ball Abutment 1mm collar	BA03
Ball Abutment 1mm collar	BA05
Ball Abutment Housing Package	BAHP1
Nylon ring Black	NRB
Nylon ring Clear	NRC
Nylon ring Pink	NRP
Spacer	SPA1
Ball Abutment Housing	BAH
Ball Abutment Analog	BAA01
Ball Abut InsertionTool	BAIT

Locator Abutment



Locator Abutment

- Simple to use low-profile solution for secure over denture attachment
- Simple nylon “male” attachments color coded for retention, simple maintenance and self aligning
- Allows for diverging angulations up to 40 degrees
- Provided in multiple tissue height orientations designed to facilitate removable implant retained overdenture



Pink (soft) 0,9 kg



White (Standard)1,5 kg




Violet (Strong)1,8 kg

Locator Abutment 0mm collar	LA00
Locator Abutment 2mm collar	LA02
Locator Abutment 4mm collar	LA04
Locator Abutment 6mm collar	LA06
Locator Abutment Housing Package	LAHP1
Locator Abutment Analog	LAA01


Full Arch Solutions: Multi-unit Abutment



St.Multi Unit




17°



30°

Flowchart – Multi unit



i-Fix-k®
Engineered for Durability



KONNECTION
|ENGINEERED TO SIMPLIFY

A New Solution for Clinicians Simplified Design & Connection for long term success



Double start threads, where thread advances the implant 1.6mm in 3.6 diameter implant & 2 mm with all other diameters with each rotation.

In contrast to i-Fix implants advance approximately 0.9mm with each rotation.

The coronal region of the i-Fix -k 4/4.5/5/6mm diameter implants are back tapered, a design that allows for maximum bone volume around the implant platform.

Implant body design The expanding tapered body is designed to attain high stability even in compromised bone situations.

Apical drilling blades facilitate insertion into minimally prepared sites.




i-Fix-k : Konnection Engineered to Simplify

	7	8	10	12	14
3.6		✓	✓	✓	✓
4.0	✓	✓	✓	✓	✓
4.5	✓	✓	✓	✓	✓
5.0	✓	✓	✓	✓	
6.0	✓	✓	✓	✓	

i-Fix-k : Master Surgical Kit

i-Fix-k
Engineered for Durability



Prefabricated Abutment

Single Unit Restoration

Screw Retained



St.MUA

Cement Retained



Cemented Abutment



CCM Hex

UCLA Hex



Angled Abutment

TiBase Hex



Screw Retained



St.MUA



CCM Non Hex

UCLA Non Hex



Multiple Unit

Cement Retained



Angled Abutment
Non Hex



Tibase Non Hex

Edentulous Solution

Screw Retained



Multiunit Abutment

Removable



Ball Abutment



Locator
Abutment

Prosthetic Components

Cover Screw

KCS01

Implant Analog

Implement the abutment on the working model.
Kind of connection - HEX 2.5



Height	Diameter	Item Code
12 MM	4.5 MM	KIA01

Healing Abutment

Diameter	Collar Height	Article Code
4.0	1	KHA4010
4.0	2	KHA4020
4.0	3	KHA4030
4.0	4	KHA4040
4.5	1	KHA4510
4.5	2	KHA4520
4.5	3	KHA4530
4.5	4	KHA4540
5.5	1	KHA5510
5.5	2	KHA5520
5.5	3	KHA5530
5.5	4	KHA5540



Prosthetic Components

Temporary Abutment

Temporary Abutment Hex KTAH

Temporary Abutment Non-Hex KTANH



Open Tray Impression Coping

For Transferring the implant position to the master model.
Available for Pick-Up(open tray) technique.

Open Tray Impression Coping Short

KOTIS

Open Tray Impression Coping Long

KOTIL



Closed Tray Impression Coping:

For Transferring the implant position to the master model.
Available for Pick-Up(close tray) technique

Closed Tray Impression Coping Short

KCTIS

Closed Tray Impression Coping Long

KCTIL



Prefabricated Final Abutments

Cement Retained Solutions



Cemented Abutment

Available in a selection of diameters and heights.
 Kind of connection -HEX 2.5 ·
 Height - 7.0mm
 Used for making the general cement type of the prosthetics.
 Recommended Tightening Torque : 30Ncm
 Use a 1.2 Hex driver

Diameter	Collar Height	Article Code
4.5	1	KCA4510
4.5	1.5	KCA4515
4.5	2.5	KCA4525
4.5	3.5	KCA4535
5.5	1	KCA5510
5.5	1.5	KCA5515
5.5	2.5	KCA5525
5.5	3.5	KCA5535

Angled Abutment

Available in a selection of diameters and heights.
 Kind of connection - HEX 2.5
 2 kinds of angle - 15°, 25° ·
 Used for adjusting the path of the prosthetics
 Use a 1.2 Hex driver.
 25° G/H 1.0 screw may protrude slightly

	Collar Height	1	2	3
15°	Hex	K15AH5001	K15AH5002	K15AH5003
	Non Hex	X	K15ANH5002	K15ANH5003

	Collar Height	1	2	3
25°	Hex	K25AH5001	K25AH5002	K25AH5003
	Non Hex	X	X	X



Prefabricated Final Abutments

UCLA Abutment

The path of abutment can be controlled through the milling process on the patient oral condition.

Material - Red (Hex) / White (Non-Hex)

Use a 1.2 Hex driver

UCLA Hex

KUCH

UCLA Non Hex

KUCNH





CCM Abutment:

Available in a selection of diameters

Material - CCM (Body) Red (Hex) White (Non-Hex)

Kind of connection - HEX 2.5

Use a 1.2 Hex driver
















Description	Article Code	
CCM Hex	KCCMH	
CCM Non Hex	KCCMNH	

Ti Base Abutment

Description	Article Code	
Ti Base abutment Hex	KTBH	
Ti Base abutment Non Hex	KTBNH	

Prefabricated Final Abutments

Multiunit Abutment

Description	Article Code	
Straight Multi Unit abut 1.5mm	KMUA15	
Straight Multi Unit abut 3.5mm	KMAU35	
17 Deg Multi Unit abut 2.5mm	K17MUA25	
17 Deg Multi Unit abut 3.5mm	K17MUA35	
30 Deg Multi Unit abut 3.5mm	K30MUA35	
30 Deg Multi Unit abut 4.5mm	K30MUA45	
Healing cap Multi Unit	HCMUA	
Abutment Screw Multi Unit	ASMUA	
Prosthetic Screw Multi Unit	PSMUA	
Temporary Titanium Coping Multi Unit	TCMUA	
Open Tray Imp Coping Multi Unit	OTICMUA	
Closed Tray Imp Coping Multi Unit	CTICMUA	
Abutment Analog Multi Unit	AAMUA	
Castable Sleeves Multi Unit Plastic	CSPMUA	
Protection Analog Multi Unit	PAMUA	

Over Denture Solutions

Ball Abutment

Over Denture solution
Use Ball abutment Driver
Maximum Torque 30Ncm.



Description	Article Code
Ball abutment 1mm	KBA01
Ball abutment 3mm	KBA03
Ball abutment 5mm	KBA05
Ball abutment Housing Package	BAHP1
Ball abutment Analog	BAA1
Ball abutment Insertion tool	BAIT01

Locator Abutment

Over Denture solution
Maximum Torque 30Ncm.



Description	Article Code
Locator abutment 0mm	KLA00
Locator abutment 2mm	KLA02
Locator abutment 4mm	KLA04
Locator abutment 6mm	KLA06
Locator abutment Housing Package	LAHP1
Locator abutment Analog	LAA01

BONE GRAFT

BCP



- 25% β tricalcium phosphate + 75% Hydroxyapatite
- Porosity 60%
- Pore size 300-500 μ m
- Mechanical resistance- 3 Mpa
- Resorption- 20% in 1.5yrs

orthopedic

dental



adbone®

BCP

75% of hydroxyapatite and
25% of beta-tricalcium phosphate



Available geometries

Granules



0,1 – 0,5mm
0,5 – 1mm
1 – 2mm
3 – 4 mm

Crunch



4 x 7mm

Cylinder



8 x 20mm

Wedges





6 x 25 x 30mm
8 x 25 x 30mm
10 x 25 x 30 mm
12 x 25 x 30mm
14 x 25 x 30mm

Blocks



4 x 4 x 20mm
5 x 10 x 15mm
8 x 8 x 20mm
15 x 15 x 20mm
15 x 20 x 30mm



New vital bone surrounding the granule of adbone®BCP, observable after 12 weeks.

High osteoclastic activity, indicating bone remodeling

TCP



- 100% β tricalcium phosphate
- Porosity 60%
- Pore size 300-500 μ m
- Mechanical resistance- 3 Mpa
- Resorption- 20% in 5-6mth

orthopedic

dental



adbone®
TCP

100% of beta-tricalcium phosphate



Available geometries

Granules



0,1 – 0,5mm
0,5 – 1mm
1 – 2mm
3 – 4 mm

Crunch



4 x 7mm

Cylinder



8 x 20mm

Wedges





6 x 25 x 30mm
8 x 25 x 30mm
10 x 25 x 30 mm
12 x 25 x 30mm
14 x 25 x 30mm

Blocks



4 x 4 x 20mm
5 x 10 x 15mm
8 x 8 x 20mm
15 x 15 x 20mm
15 x 20 x 30mm

12 week follow-up



Due to our interconnected structure and large pore size, we can see new bone formation in the middle of the TCP particule

BCX



INDICATIONS

- Full mouth rehabilitation
- Graftless solution
- Thin ridge situation
- Bicortical engagement

	A	B	C	D	E	F	G	H
3.6x8	1	4.5	3.5	1.5	3.1			1.6
3.6x10	1	6	4	1.5	3.1			1.6
3.6x12	1	6	6	1.5	3.1			1.6
3.6x14	1	9	5	1.5	3.1			1.6
3.6x17	1	9	8	1.5	3.1	3	2.5	1.6
3.6x20	1	9	11	1.5	3.1	3	2.5	1.6
3.6x23	1	9	14	1.5	3.1	3	2.5	1.6
3.6x26	1	9	17	1.5	3.1	3	2.5	1.6

A	—	Tip Length
B	—	Thread Length
C	—	Shaft Length
D	—	Apical Core Diameter
E	—	Apical Major Diameter
F	—	Collar Length
G	—	Collar Diameter
H	—	Pitch

PROSTHETIC OPTIONS



Impression Coping



BCX Analog



BurnOut Cap



INDICATIONS

- Immediate prosthetic loading
- Elimination of screw loosening
- Reduces cost & time

	A	B	C	D	E	F	G	H	I	J
3510	0.65	10	10	2	1.8	0.9	3.5	2.46	1.55	0.5
3512	0.65	12	8	2	1.8	0.8	3.5	2.46	1.55	0.5
3514	0.65	14	7	2	1.8	0.8	3.5	2.46	1.55	0.5
3516	0.65	16	6	2	1.8	0.8	3.5	2.46	1.55	0.5
4008	0.7	8	13.5	2	2.2	1.2	4	2.95	1.55	0.5
4010	0.7	10	13	2	1.8	0.8	4	2.95	1.55	0.5
4012	0.7	12	11	2	1.8	0.8	4	2.95	1.55	0.5
4014	0.7	14	9	2	1.8	0.8	4	2.95	1.55	0.5
4016	0.7	16	8	2	1.8	0.8	4	2.95	1.55	0.5
4508	0.75	8	12.5	2.3	2.8	1.4	4.5	3.05	1.55	0.7
4510	0.75	10	13	2.3	2.3	0.9	4.5	3.05	1.55	0.7
4512	0.75	12	10.5	2.3	2.3	0.9	4.5	3.05	1.55	0.7
4514	0.75	14	9	2.3	2.3	0.9	4.5	3.05	1.55	0.7
5008	0.85	8	14	2.3	3.2	1.8	5	3.55	1.55	0.7
5010	0.85	10	14.5	2.3	2.6	1.2	5	3.55	1.55	0.7
5012	0.85	12	12	2.3	2.6	1.2	5	3.55	1.55	0.7
5014	0.85	14	11	2.3	2.3	0.9	5	3.55	1.55	0.7
5508	0.95	8	18	2.5	3.2	1.8	5.5	4.04	1.55	0.7
5510	0.95	10	18	2.5	2.5	1.1	5.5	4.04	1.55	0.7
5512	0.95	12	15	2.5	2.5	1.1	5.5	4.04	1.55	0.7

A	—	Tip Length
B	—	Thread Length
C	—	Total Angle
D	—	Shaft Diameter
E	—	Apical Diameter
F	—	Apical Core Diameter
G	—	Implant Crestal Diameter
H	—	Implant Crestal Core Diameter
I	—	Pitch
J	—	Thread Depth

PROSTHETIC OPTIONS



Impression Coping



CIX Analog



BurnOut Cap

Thank You

